

digital hotplates,
stirrers and
hotplate stirrers
are the ideal
solution for
demanding users
in all laboratory
environments



CLEAVER SCIENTIFIC LTD

Unit 41, Somers Road Industrial Estate,
Rugby, CV22 7DH
United Kingdom

T_ +44 (0)1788 565300

E_ INFO@CLEAVERSCIENTIFIC.COM

W_ WWW.CLEAVERSCIENTIFIC.COM

Digital Hotplates and Stirrers

- Large, backlit LCD display
- Ceramic work surface, 6.5 x 6.5 in.
- Safety LED indicates hot surface
- Control actual temperature (with optional probe)
- Three models: heat-stir, heat only or stir only

With a durable and chemically resistant ceramic surface, Cleaver Scientific's digital hotplate, stirrer and hotplate stirrer are the ideal solution for demanding users in all laboratory environments. The minimal footprint (18 x 26 cm) allows for use in crowded spaces such as fume hoods while the 16.5 cm square plate makes these units compatible with a wide range of commonly used vessels such as beakers, bottle and conical flasks.

Fast and precise adjustment of speed and temperature is achieved with advanced microprocessor technology, and a large backlit LCD display offers easy viewing of current parameters. A safety LED indicates temperatures over 50°C. An optional external thermometer and support rod allows temperature control of the sample by direct feedback to the microprocessor, maintaining temperature to within $\pm 0.5^{\circ}\text{C}$.



TECHNICAL SPECIFICATIONS

TEMP. RANGE:	AMBIENT +5 TO 380°C
SPEED RANGE:	200 TO 1500 RPM
PLATFORM DIMENSIONS:	6.5 x 6.5 in/ 16.6 x 16.5 cm
OPERATING TEMP. RANGE	+4 TO +65°C
DIMENSIONS: (W X D X H)	7 x 10.5 x 4 in. 18 x 26 x 10.1 cm
WEIGHT:	5 LBS (2.2kg)
ELECTRICAL:	120 OR 230V, 50-60Hz 500W
WARRANTY:	2 YEARS



ORDERING INFORMATION

Cat. No.	Description	Accessories	
CSL-DHOTPLATE	Digital Hotplate, 6.5 x 6.5 in. - 230V	TEMPROBE	Optional Temperature probe.
CSL-DHOTPLATE\$	Digital Hotplate, 6.5 x 6.5 in. - 110V	SUPPROD	Optional Support Rod
CSL-DSTIR	Digital Magnetic Stirrer, 6.5 x 6.5 in.- 230V		
CSL-DSTIR\$	Digital Magnetic Stirrer, 6.5 x 6.5 in.- 110V		
CSL-DHOTSTIR	Digital Hotplate Magnetic Stirrer 6.5 x 6.5 inch. - 230V		
CSL-DHOTSTIR\$	Digital Hotplate Magnetic Stirrer 6.5 x 6.5 inch. - 110V		

CLEAVER SCIENTIFIC LTD

Unit 41, Somers Road Industrial Estate,
Rugby, CV22 7DH
United Kingdom

T_ +44 (0)1788 565300

E_ INFO@CLEAVERSCIENTIFIC.COM

W_ WWW.CLEAVERSCIENTIFIC.COM